

Newborn Screening and Diagnosis of Fatty Acid Oxidation Disorders

Carla Cuthbert, PhD, FACMG, FCCMG

Chief, Newborn Screening and Molecular Biology Branch
National Center for Environmental Health
Centers for Disease Control and Prevention

Newborn Screening

A Public Health Program

- Aimed at identification of conditions for which early intervention can prevent mortality, morbidity, and disabilities
- Performed by analysis of diagnostic markers in blood spots collected on filter paper at birth
- Largest population-based genetic screening effort in the US: ~ 4 Million babies screened annually

The Newborn Screening System

- **Screening**
 - Laboratory analysis of newborn bloodspot
- **Follow-up of an abnormal result**
 - Rapid location, follow-up and referral of the screen-positive infant
- **Diagnosis**
 - Clinical and biochemical evaluation to diagnose or rule out the disorder
- **Management**
 - Rapid planning and implementation of long-term therapy
- **Evaluation**
 - Assessment of the NBS program: efficiency of follow-up & intervention, benefit to society



**Alaska
Newborn Screening
Program**

Newborn Screening: Short Term Follow-up

Abnormal result obtained from the State Lab



CMS nurse from follow-up team notifies Genetics on-call Physician about the Abnormal Newborn Screen result



Genetics on-call Physician locates patient, makes clinical assessment of patient status – recommends clinic visit



**Infant sample collected (Plasma, Urine, Whole Blood)
Diagnostic evaluation performed**



Positive cases: long-term patient management & ongoing care

Follow-up Diagnostic Testing

- **Routine Chemistries**

- Electrolytes, blood gases, anion GAP
- Blood sugar, ketones, lactate, ammonia

- **Biochemical Genetics Testing**

- Plasma acylcarnitine (Green top tube – 0.5 ml min whole blood)
- Carnitine Status (Red top tube – 0.5 ml min whole blood)
- Urine organic acids (No preservatives – 3 ml min vol)
- Plasma amino acids (Green top tube – 0.5 ml min whole blood)

- **Other specialized testing**

- Enzyme analysis (eg Fatty acid oxidation Probe Assay)
- Mutation testing

Tandem Mass Spectrometry

Clinical Chemistry 49:11
1797–1817 (2003)

Review

Use of Tandem Mass Spectrometry for Multianalyte Screening of Dried Blood Specimens from Newborns

DONALD H. CHACE,* THEODORE A. KALAS, and EDWIN W. NAYLOR

***Tandem Mass
Spectrometer***



***A highly sensitive and
specific tool for diagnostic
sample analysis***


Tandem Mass Spectrometry

- Identification of unknown compounds
- Quantitation of known analytes using stable isotopes
- Provides structural and chemical information of molecules
- High sensitivity and specificity
- Fast analytical time
- Potential for automation for high throughput tests

(modified from ASMS, 1989)

Acylcarnitine Analysis by MS/MS

An Acylcarnitine molecule has the following structure:



Members of this family of compounds have varying lengths of their R-Group

C3 acylcarnitine ... 3 carbon units
C8 acylcarnitine ... 8 carbon units
C16 acylcarnitine ... 16 carbon units

- Specific Acylcarnitines accumulate in fatty acid oxidation disorders (and certain organic acidemias) and form a distinct “pattern”
- Acylcarnitine analysis is very important in: prenatal diagnosis, newborn screening, evaluation of symptomatic patients, and postmortem screening

Newborn Screening by MS/MS

Disorders of fatty acid oxidation

2,4-Dienoyl-CoA reductase deficiency
Carnitine acylcarnitine translocase deficiency
Carnitine palmitoyltransferase I deficiency
Carnitine palmitoyltransferase II deficiency
Carnitine transport defect
Electron transfer flavoprotein deficiency
ETF ubiquinone oxidoreductase deficiency
Long-chain L-3-OH acyl-CoA dehydrogenase def.
Medium-chain acyl-CoA dehydrogenase deficiency
Medium-chain L-3-OH acyl-CoA dehydrogenase def.
Medium chain ketoacyl-CoA thiolase deficiency
Short-chain acyl-CoA dehydrogenase deficiency
Trifunctional protein deficiency
Very long-chain acyl-CoA dehydrogenase def.

Disorders of amino acid metabolism

Arginase deficiency
Argininosuccinate lyase deficiency
Argininosuccinate synthase deficiency
Maple syrup urine disease (MSUD)
Citrin deficiency
Cystathionine β -synthase deficiency
Methionine adenosyltransferase deficiency
Mitochondrial ornithine transport defect (HHH)
Phenylalanine hydroxylase deficiency (PKU)
Defects of bipterin metabolism
Fumarylacetoacetase deficiency
Tyrosine aminotransferase deficiency

Collective incidence


1:2-4,000 newborns

**Clinical impact (US):
Up to 2,000 cases/year**


Disorders of organic acid metabolism

2-CH₃ butyryl-CoA dehydrogenase deficiency
2-CH₃ 3-OH butyryl-CoA dehydrogenase deficiency
3-OH 3-CH₃ glutaryl-CoA lyase deficiency
3-CH₃ crotonyl-CoA carboxylase deficiency
3-CH₃ glutaconyl-CoA hydratase deficiency
Isobutyryl-CoA dehydrogenase deficiency
Isovaleryl-CoA dehydrogenase deficiency
Glutaryl-CoA dehydrogenase deficiency
Malonyl-CoA carboxylase deficiency
Methylmalonyl-CoA mutase deficiency
Disorders of cobalamin metabolism
 β -ketothiolase deficiency
Multiple carboxylase deficiency
Propionyl-CoA carboxylase deficiency

Normal Acylcarnitine Profile




VLCAD Acylcarnitine Profile



Normal


**Free
Carnitine**




**Patient
(on MCT diet)**

**Normal Profile in known
VLCAD patient on diet**


**Free
Carnitine**




Postmortem Bile (LCHAD)



Postmortem DBS (LCHAD)




NBS DBS (LCHAD)



Fatty Acid Oxidation Probe Assay

- **Monolayer of skin fibroblasts / amniocytes**
- **Growth media removed and replaced with “*in vitro* probe” media containing palmitic acid and L-carnitine**
- **Incubated for 72 hours, media is collected and subjected to acylcarnitine analysis**
- **Adherent cells are harvested - total protein analysis**




Enzyme Assay for VLCAD - Results

	Control Amniocytes	VLCAD Patient Amniocytes
MCAD Activity (mU/mg Protein)	1.01	1.21
VLCAD Activity (mU/mg Protein)	1.52	Not detected

	Control Fibroblasts	VLCAD Patient Fibroblasts
MCAD Activity (mU/mg Protein)	1.31	1.07
VLCAD Activity (mU/mg Protein)	3.61	Not detected

Urine Organic Acids

- **Organic Acids – What are they?**
 - Water-soluble compounds containing one or more carboxylic group as well as other functional groups (-keto, -hydroxy)
 - Intermediate metabolites of all major groups of organic cellular components: amino acids, lipids, carbohydrates, nucleic acids and steroids
- **Accumulation of metabolites which are not present under physiological conditions**




Urine Organic Acids

- **Urine Sample Collection**
 - Urine should be collected without preservatives
 - Collection should avoid fecal contamination
 - Frozen immediately
- **Sample Preparation and Analysis**
 - Liquid-liquid extraction of acidified urine into organic solvent
 - Evaporation of solvent and derivatization of residue
 - Gas chromatography – Mass spectrometry analysis
- **Specific diseases have characteristic organic acid elevations**

Organic Acid Analysis

Propionic Acidemia



Plasma Amino Acids

- **Plasma Sample Collection**
 - Timely centrifugation and separation of plasma specimens is critical to prevent artifacts
 - Must be refrigerated for the short term (< 4 hrs.) or frozen (-20°C) to arrest amino acid deterioration
 - Hemolyzed samples are not viable
- **Sample Preparation and Analysis**
 - Protein precipitation and centrifugation
 - Ion exchange chromatography to separate
 - Post-column derivatization & detection at 440 to 570 nm
- **Specific diseases have characteristic amino acid elevations**

Plasma Amino Acids Analysis

Tyrosinemia

